

Gloway Invent Oy

Dental Vacuum autoclave

MODEL: BTD17L/BTD23L



G

Gloway Invent Oy

Product detail

Functional Character

- · Advanced steaming pressure way of sterilization recognized internationally.
- ·Operating menus convenience, consumer may setup day, time etc.
- •The structure of CROSS generational full automatic electronic lock.
- ·Deep-stretched shape high quality stainless steel internal bladder
- ·High-luminance liquid crystal display
- ·Built-in temperature sensing printing, and equipped with pc interface, operating with remote maintenance.
- ·Setup warm-up key, in order to keep room of sterilization temperature around 70°C .
- ·Built-in and set by itself procedure of antisepsis, and the range of the antisepsis is wide
- ·Built-in cleaning and sewage equipment
- ·When outlaying it can connect with the system of distilled water, it can also connect with the sewage pipe directly.

Application Range

Puh: 050 323 9595

Pressurized steam sterilizer is mainly applied to operating room. Stomatological Department, Ophthalmologic Department and biomedicine research unit, and it can sterilize some substances such as instruments of operation, instruments of Stomatological Department, dressing, glassware, injecting apparatus.